

FOR *in vitro* USE ONLY

Reconstituted

in vitro translation kit

PURESYSTMTEMTM originator

Instruction Manual

originator8	#Por008
originator96	#Por096
originatorS-S8	#PoS008
originatorS-S96	#PoS096

Postgenome Co., Ltd.

9-1 Azanaka Kiriura, Nishiitsushiro, Ichinomiya, Aichi JAPAN
Tel: +81-586-59-4385 Fax: +81-586-58-2339

IMPORTANT NOTICE to Purchaser

- This kit cannot be used in human subjects, in clinical trials or for diagnostic purposes involving human subjects.
- Contents of this kit are subject to change without notice.

Safety Instructions

Please note that Sol. A and Sol. B contains <1% DTT (Dithiothreitol). No other hazardous substances are included in reportable quantities. The synthesized protein may have strong toxicity or show pathogenic properties. To insure your safety, please take all the necessary precautions such as using the safety cabinet or protecting your eyes and skin by use of eyewear and gloves at all times. If your skin comes in direct contact with any of the reagents, or if you have swallowed any reagents by mistake, please obtain appropriate medical assistance. Please also insure to inactivate all unused reagents by autoclave or dry-heat sterilization before disposal. If your institution or country holds its own laws or regulations concerning experimental procedures, material transfer, export and disposal of the materials, please follow the instructions as stated.

Table of Contents

1. Quick Protocol for Protein Synthesis Reaction
2. Kit Components
3. Storage Conditions
4. Description
 - 4-1 Background of PURESYSTEM
 - 4-2 PURESYSTEM originator
5. Procedures
 - 5-1 Preparation of Template DNA
 - 5-2 Protein Synthesis Reaction
 - 5-3 Purification of Synthesized Protein-Quick Purification Protocol
 - 5-4 Analysis of Synthesized Protein
6. Synthesis and Quick Purification of DHFR(Positive Control)
7. Trouble shooting
8. Technical Service
9. Related Products

1. Quick Protocol for Protein Synthesis Reaction

1. Thaw Sol. A and Sol. B on ice.
2. Assemble in new reaction tube(s) as follows;

<input type="checkbox"/>	Nuclease-free water	μl
<input type="checkbox"/>	Sol. A	25 μl
<input type="checkbox"/>	Sol. B	10 μl
<input type="checkbox"/>	Template DNA	μl^{*1}
<hr/>		
	Total	50 μl

*1; 1 pmol PCR product or 0.5 μg plasmid DNA

3. Mix gently and centrifuge briefly.
4. Incubate the reaction(s) at 37°C for 1 hour.
5. Stop the reaction by placing the tube(s) on ice.
6. Analyze the synthesized protein(s).

2.Kit Components

Please check that your kit contains the reagents listed below.

2-1. PURESYSTMTEM originator 8(prod.#Por008)

Each kit contains sufficient reagents for 8 x 50 µl reactions.

Label (Tube Color)	Amount
Sol. A (Yellow)	8x 25µl
Sol. B (Red)	8x 10µl
Temp. DHFR (Blue)	1 x 5µl
Univ. Primer (Clear)	1 x 50µl

2-2. PURESYSTMTEM originator 96(prod.#Por096)

Each kit contains sufficient reagents for 96 x 50 µl reactions.

Label (Tube Color)	Amount
Sol. A (Yellow)	96 x 25 µl
Sol. B (Red)	96 x 10 µl
Temp. DHFR (Blue)	2 x 5 µl
Univ. Primer (Clear)	2 x 500 µl

2-3. PURESYSTMTEM originator S-S8(prod.#PoS008)

Each kit contains sufficient reagents for 38x 50 µl reactions.

Label (Tube Color)	Amount
Sol. A (Yellow)	8x 25µl
Sol. B (Red)	8x 10µl
Temp. DHFR (Blue)	1 x 5µl
Univ. Primer (Clear)	1 x 50µl

2-4. PURESYSTMTEM originator S-S96(prod.#PoS096)

Each kit contains sufficient reagents for 96 x 50 µl reactions.

Label (Tube Color)	Amount
Sol. A (Yellow)	96 x 25 µl
Sol. B (Red)	96 x 10 µl
Temp. DHFR (Blue)	2 x 5 µl
Univ. Primer (Clear)	2 x 500 µl

- Sol. A and Sol. B must be thawed just before use on ice and avoided multiple freeze-thaw cycles.
- Temp. DHFR is plasmid DNA containing dihydrofolate reductase (DHFR) gene (0.2 µg/µl) and can be used as positive control. It is added directly to the synthesis reaction.
- Univ. Primer is a forward primer for second-step PCR and its concentration is 2 µM.

3.Storage Condition

Store all the components in the kit at -80°C. Avoid multiple freeze-thaw cycles.

4.Description

4-1. Background of PURESYSTEM Technology

PURESYSTEM is a novel reconstituted *in vitro* translation system which consists of about 30 purified enzymes necessary for transcription, translation and energy recycling.* In **PURESYSTEM**, all the factors for transcription and translation are tagged with a hexahistidine; 3 initiation factors (IF1, IF2, IF3), 3 elongation factors (EF-G, EF-Tu, EF-Ts), 3 release factors (RF1, RF2, RF3), ribosome recycling factor, 20 aminoacyl-tRNA synthetases, methionyl-tRNA formyltransferase and T7 RNA polymerase. The reagents also contains *E. coli* 70S ribosome, amino acids, NTPs, *E. coli* tRNA and energy recycling system, so the target protein can be synthesized just by addition of the template DNA to the reaction.

* Shimizu Y., *et al.*, (2001) *Nature Biotechnology*, vol.19, p.751-755.

4-2. PURESYSTEM originator

PURESYSTEM originator is a kit based on PURESYSTEM technology

High Purify

Because **PURESYSTEM** is a reconstituted system in which all the components are identified, the reagents in **PURESYSTEM** originator are essentially free from unidentified components. Therefore, the synthesized product is hardly susceptible to unwanted degradation or modification.

No Tags Required

Because all the proteinous factors are His-tagged and no unknown components are present within the reaction, the synthesized protein can be easily purified by removal of the His-tagged factors and ribosome using metal affinity resin and ultrafiltration membrane, respectively (Figure 1). Therefore, no unexpected effects of tags needs to be considered (e.g.; misfolding or aggregation), promising a hassle-free purification of the synthesized protein.

Time Saving

All you need is to follow the 3 steps as shown in Figure 1. It takes only 10 minutes for your handling and 3 hours for incubation and centrifugation.

4-3. Related Kits

PURESYSTEM™ custom

PURESYSTEM is a reconstituted system, so the reagent composition can be easily modified according to the purpose of experiment. This can be ordered as **PURESYSTEM custom**, which is made upon request.

PURESYSTEM™ S-S

PURESYSTEM S-S is the *in vitro* transcription/translation system which allows the synthesized protein to form correct disulfide bond(s). The high purity and flexibility of **PURESYSTEM** enables the redox conditions of reaction mixture to be reproducible.




Figure 1.

Figure1.

Schematic diagram of protein systhesis and purification
by PURESYSYSTEM

5.Procedures

5-1. Preparation of Template DNA

5-1-1. Considerations for Template DNA Construction

Both PCR product and plasmid DNA can be used as the template DNA for protein synthesis reaction in **PURESYSYSTEM**. The template DNA needs to contain the following;

- start codon (ATG)
- stop codon (TAG, TGA or TAA)
- T7 promoter sequence at upstream from the coding sequence
- Ribosomal binding site (RBS; SD sequence) at approximately 10 bp upstream from start codon
- Additional sequence of 6 bp or more at downstream from stop codon (for PCR product)
- Terminator sequence at downstream from stop codon (for plasmid DNA)

We recommend PCR product as the template DNA due to the following reasons;

- PCR product can be used without further purification.
- No need for time-consuming construction of plasmid DNA suitable for protein synthesis.
- Multiple templates can be prepared at a time.
- The coding sequence can be easily modified.
eg) deletion, point mutation, division into several domains and addition of tag sequence.

5-1-2. Generation of the Template DNA by 2-step PCR

The template DNA can be generated by 2-step PCR (Figure 2). In the first (1st) PCR, the adapter sequence is added at 5'-terminus of the gene of your interest, and all the regulatory sequences are added in the second (2nd) PCR. It is possible to amplify the template DNA by 1-step PCR using forward primer containing 5'-terminal region of the coding sequence and all regulatory sequences.

1. Design the forward and reverse primers as shown in Table 1.

2. Perform 1st PCR from DNA containing the gene of interest using these primers and high-fidelity polymerase.
3. Analyze the 1st PCR product by agarose gel electrophoresis. When detected as a single band, the 1st PCR reaction mix can be used in the second step without further purification. If any additional bands are detected, please optimize the PCR conditions or excise the band of the product from agarose gel.
4. Perform 2nd PCR from the 1st PCR product using universal primer included in the kit and reverse primer used at the 1st PCR.
5. Analyze the 2nd PCR product by agarose gel electrophoresis. The 2nd PCR product will be approximately 70 bp longer than the 1st PCR product, and this can be used for protein synthesis reaction without further purification.

Table1. Primer List

primer	sequence
Forward primer	5'-AAGGAGATATACCA-ATG-N ₁₄₋₂₀ -3'
Reverse primer	5'-TATTCA-TTA-N ₁₄₋₂₀ -3'
Universal primer (Univ. Primer)	5'-GAAATTAATACGACTCACTATAGGGAGACCACAACGGTTTCC CTCTAGAAATAATTTTGTTTAACTTTAAGAAGGAGATATACCA-3'

5-1-3. Considerations for plasmid DNA


Plasmid DNA can also be used as the template DNA for **PURESYSYSTEM**. The plasmid DNA will need to fulfill the following elements and structural features;

- Any commercially available expression vectors can be used which contains T7 promoter sequence, RBS (approximately 10 bp upstream from the initiation codon) and terminator sequence.
- Some commercially available plasmid purification kits are not suitable for preparation of template DNA for **PURESYSYSTEM**. Please contact technical service for details.
- The plasmid DNA needs to be dissolved in RNase-free water or buffer without EDTA, in order to avoid inhibition of protein synthesis reaction.
- RNA needs to be removed. When RNase is used for the removal of RNA, complete inactivation of RNase is necessary by Phenol treatment or by similar inactivation

treatments.

- Quick purification method described in this manual cannot be used when the target gene is introduced to expression vector with a histidine-tag sequence.

1st Step PCR



2nd Step PCR




Figure2. Generatrion of Template DNA by 2-step PCR

5-2. Protein Synthesis Reaction

5-2-1. Precautions concerning reagents and experimental procedures

If you are unfamiliar with an *in vitro* translation, we recommend you to check the experimental procedures by using positive control (DHFR) included in the kit.

In order to avoid nuclease contamination, please use Nuclease-free tubes and tips, and wear gloves and mask.

5-2-2. Procedures

1. Thaw the necessary number of Sol. A and Sol. B on ice.

Please avoid multiple freeze-thaw cycles, because some enzymes in the reagents may be inactivated.

2. Assemble in new reaction tube(s) as follows;

<input type="checkbox"/> Nuclease-free water	μl
<input type="checkbox"/> Sol. A	25 μl
<input type="checkbox"/> Sol. B	10 μl
<input type="checkbox"/> Template DNA	μl* ¹
<hr/>	
Total	50 μl

*1; 1 pmol PCR product or 0.5 μg plasmid DNA

When labeled amino acid is incorporated into synthesized protein, add the labeled amino acid just before protein synthesis. For example, add 1 μl of 0.4 MBq/μl [³⁵S]methionine before addition of template DNA.

3. Mix gently and centrifuge briefly.
4. Incubate the reaction(s) at 37°C for 1 hour.
5. Stop the reaction by placing the tube(s) on ice.

The reaction mixture is now ready to be analyzed. (e.g., purification, electrophoresis, enzyme assay, binding assay).

Some proteins may require couple of hours for correct folding to take place after synthesis reaction. If correct folding is known to be crucial for your protein of interest, please place the tube in the incubator.

Refer to section 5-3 for purification.

If purification step is unnecessary, the reaction mixture may be directly used for your experiment.

5-3. Purification of Synthesized Protein - Quick Purification Protocol

If your proteins are characterized and you are familiar with the purification methods for your proteins, please proceed with the methods of your choice.

If the purification methods of the synthesized protein are unknown and the molecular weight is less than 100 kDa, the protein may be purified by the quick purification method as follows.

5-3-1. Additional Equipments necessary for Quick Purification

Ni-NTA Agarose (QIAGEN)

Nanosep 100K (Pall)

Micro Bio-Spin Chromatography Column (BIO-RAD)

Vortex Mixer

5-3-2. Procedure 1

1. Add equal volume of H₂O to the reaction.
2. Apply the diluted reaction mixture to Nanosep 100K.
3. Centrifuge at 1,500xg at 4°C for 30 minutes.
4. Transfer the permeate to a new tube (preferably 2 ml round-bottom tube).
5. Add 0.1 vol. of Ni-NTA Agarose.
6. Vortex for 1 hour at 4°C vigorously.
7. Apply the mixture to the Micro Bio-Spin Chromatography Column.
8. Centrifuge at 1,500xg at 4°C for 1 minute.
9. Recover the permeate.

5-3-3. Procedure 2

This procedure is recommended when the synthesized protein could be purified according to the Procedure 1.

1. Add equal volume of H₂O to the reaction.
2. Transfer the diluted reaction mixture to a new tube (preferably 2 ml round-bottom tube).
3. Add 0.1 vol. of Ni-NTA Agarose.

4. Vortex for 1 hour at 4°C vigorously.
5. Apply the mixture to Nanosep 100K.
6. Centrifuge at 1,500xg at 4°C for 30 minutes.
7. Recover the permeate fraction.

5-4. Analysis of Synthesized protein

Synthesized proteins can be detected and analyzed by methods such as SDS-PAGE followed by gel staining or western blotting and enzyme activity assay. In this section, a general method for SDS-PAGE and gel staining is shown.

If protein synthesis reaction is carried out in the presence of radiolabeled amino acids, the product can be analyzed by fluorography.

1. Mix 1-10 μ l of samples with SDS-PAGE sample buffer.
2. Load the samples to SDS-PAGE gel and carry out electrophoresis.
3. Stain the gel by Coomassie Brilliant Blue, Silver, or Fluorescent dye such as SyproOrange.


The yield of the target product may be dependent on each protein. If the yield is low, please concentrate the samples by methods such as trichloroacetic acid and acetone. Depending on the staining method or kit, tRNA band may be detected around 20-30 kDa. The background can be reduced by treatment with RNase after synthesis reaction. For example, add 3.75 μ g of RNaseA in 50 μ l of reaction and incubate at 37°C for 15 min.

6. Synthesis and Quick Purification of DHFR (Positive Control)

1. Thaw Sol. A, Sol. B and Temp. DHFR on ice.
2. Assemble in a new reaction tube as follows;

Nuclease-free water	12.5 μ l	
Sol. A	25 μ l	
Sol. B	10 μ l	
Temp. DHFR	2.5 μ l	(0.5 μ g)
<hr/>		
Total	50 μ l	

3. Mix gently and centrifuge briefly.
4. Incubate the reaction at 37°C for 1 hour.
5. Stop the reaction by placing the tube on ice.
6. Add 50 μ l of H₂O to the reaction.
7. Apply the diluted reaction mixture to Nanosep 100K.
8. Centrifuge at 1,500xg at 4°C for 30 minutes.
9. Transfer the permeate to a new 2 ml round-bottom tube.
10. Add 10 μ l of Ni-NTA Agarose.
11. Vortex for 1 hour at 4°C vigorously.
12. Apply the mixture to the Micro Bio-Spin Chromatography Column.
13. Centrifuge at 1,500xg at 4°C for 1 minute.
14. Recover the permeate.
15. Mix 10 μ l of each sample (equiv. 5 μ l of the synthesis reaction) with SDS-PAGE sample buffer.
16. Incubate the samples at 95°C for 5 minutes.
17. Load to 12.5 % SDS-PAGE gel.
18. Carry out electrophoresis.
19. Stain the gel by Coomassie Brilliant Blue R-250.



1. MW marker
2. after Synthesis (- template DNA)
3. after Synthesis
4. Gel-filtrate
(after removal of Ribosome)
5. after Removal of Ni-NTA Agarose
(after removal of His-Tagged Factors)

Figure 3. Synthesis and quick purification of DHFR

7. Troubleshooting

7-1. Preparation of Template DNA

No product at the 1st step PCR

Primer is not optimal.

→ Redesign the primer.

Condition is not optimal.

→ Check reaction conditions such as annealing temperature, extension time and number of cycles.

→ Perform the PCR reaction in the presence of additives such as DMSO.

No product at the second-step PCR

Primer is not optimal.

→ Redesign the primer.

→ Prepare the template DNA by 1 step PCR using the long primer as below.

5'-GAAATTAATACGACTCACTATAGGGAGACCACAACGGTTTCC

CTCTAGAAATAATTTTGTTTAACTTTAAGAAGGAGATATACCA-ATG-N₁₄₋₂₀-3'

Condition is not optimal.

→ Check reaction conditions such as annealing temperature, extension time and number of cycles.

→ Perform the PCR reaction in the presence of additives such as DMSO.

7-2. Protein Synthesis by PURESYSTEM

No control is synthesized.

Some components in the kit are inactivated.

→ Store the reagents in the kit at -80 °C and avoid multiple freeze-thaw cycles.

The reaction is contaminated with nuclease.

→ Wear gloves and mask and use nuclease-free tubes and tips.

Control is synthesized, but the target is not.

DNA sequence of the template DNA is incorrect.

→ Check the sequence of the template DNA.

Purity and/or concentration of the template DNA is not optimal.

→ Check purity and concentration of the template DNA.

The reaction is contaminated with nuclease.

→ Wear gloves and mask and use nuclease-free tubes and tips.

Transcript may have formed a rigid secondary structure.

- Introduce mutations in the sequence to avoid secondary structure formation.
- Insert a tag coding sequence at the 5' terminus of the target sequence.

Products (transcripts and/or polypeptides) inhibit protein synthesis reaction.

- Please synthesize control protein (DHFR) in the same reaction mixture and check the yield of DHFR.

The target protein is synthesized, but the yield is very low.

Purity and/or concentration of the template DNA is not optimal.

- Check purity and concentration of the template DNA.

The reaction is contaminated with nuclease.

- Wear gloves and mask and use nuclease-free tubes and tips.

Transcript may form a rigid secondary structure.

- Introduce mutations in the sequence to avoid secondary structure formation.

The target protein is synthesized, but it is insoluble.

Products form aggregates.

- Protein solubility is improved when protein synthesis reaction is performed at lower temperature than 37°C.

8. Technical Service

For information about **PURESYS**TEM technology and kits (including **custom** kit)

info@post-genome.com

The use of this kit does not provide the purchaser in any way, expressed or implied, under any patents, patent applications, or other proprietary rights by BioComber Co., Ltd. In particular, no expressed or implied licenses or other rights are provided to use this kit for any commercial purposes, and the use of this kit is strictly limited to internal research activities of the purchaser. For information concerning availability of commercial license, please contact BioComber Co., Ltd.